# SINTERFILTER HIGH GLASS

**ROBU**<sup>®</sup>  
GLASFILTER-GERÄTE

2011/2012

# ROBU

## Der Name

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# VitraPOR®

## The Name

**VitraPOR®** is the trade-name for a complete series of glass filter-products, equipped with filter elements. ROBU-products represent the internationally established type of borosilicate glass 3.3 according to DIN-ISO 3585 and english, american and french standards. Highest possible chemical resistance, minimum thermal expansion as well as the high thermal shock resistance are the characteristic points of ROBU glass filter products. Due to these features the ROBU filters are specially designed for use in the chemical, biotechnical, pharmaceutical and laboratory field.

Standard fritted discs are available:

- in round shape from 5 mm to 380 mm diameter
- in round shape bi-plane ground
- in round shape bi-plane ground and with fused edge
- filter candles in cylindrical or conical shapes from 9 mm to 40 mm diameter
- all filterapparatus and absorbers

Custom made **VitraPOR®**-filters:

- any shape up to 400 mm in almost any thickness
- all custom made discs are available unground, bi-plane ground and with fused edge
- any cylindrical or conical filter candle according to your design up to 500 mm length

Please contact us for any custom made item. We are the experts for any special design and are able to solve your problems fast and at reasonable costs.

**Prices in EURO/piece ex works according to our terms of delivery and payment.  
VAT will be invoiced according to European regulations.  
All previous prices are no longer valid.**

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
## Table of contents

Article	Page
Carrier substrate (for cell cultivation)	4
Filterdiscs, centred	4, 5
Filterdiscs, bi-plane	6, 7
Discs for washing bottles	8
Filtercandles	8, 9
Vacuum-/Gas distribution filter	9
HCL Absorber CSB-Reaction vessel	9
Gas distribution tubes	10
Filter crucibles	10
Filter funnels and accessories	11
Filter funnels, conical	12
Pipeline filters	12
Immersion tubes	12
Micro immersion tubes	12
Micro filter funnels	12
Filter tubes, Allihn	13
Chromatography-Columns	13
Suction flasks	13
Gas washing bottles	14, 15
Filter crucibles, Gooch	16
Gas filter funnels, conical	16
Filter crucibles, with mantle	16
VitraPOR SPE Glass columns, empty	16
Technical Data	17, 18, 19, 20
Porosities	21
Flow rate and drop in pressure	22
Safety advice	23, 24
Specialities Twin-filter Casting Technologies Porous Glass Beads	25
Terms of delivery and payment	26

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
**Porous Glass Beads**

Substrate for Cell Cultivation  
Pure Borosilicate Glass 3.3

Best.-Nr.	€/Stück	Por.	B	VE
Cat. -No.	€/piece		mm	

80 06 8	483,15	2	4,0	1 Liter
---------	--------	---	-----	---------

80 08 0	2730,00	2	1,5	1 Liter
---------	---------	---	-----	---------

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**Filter discs**

edge not fused, centred

\*The faces of our filterdiscs from the 15 and 51 – 53 series are press moulded and may show thickness deviations. Please refer to the discs of our 16 and 61 – 63 series on the following pages for more specific requirements.


Best.-Nr.	€/Stück	Por.	B	*S	VE
Cat. -No.	€/piece		mm	ca. mm	

15 05 0	2,90	0	5	2,8	10
15 05 1	2,90	1			10
15 05 2	2,90	2			10
15 05 3	2,90	3			10
15 05 4	3,30	4			10
15 05 5	9,65	5			10

15 10 9	2,90	00	10	2,8	10
15 10 0	2,90	0			10
15 10 1	2,90	1			10
15 10 2	2,90	2			10
15 10 3	2,90	3			10
15 10 4	3,30	4			10
15 10 5	9,65	5			10

15 20 9	3,00	00	20	3,0	10
15 20 0	3,00	0			10
15 20 1	3,00	1			10
15 20 2	3,00	2			10
15 20 3	3,00	3			10
15 20 4	3,35	4			10
15 20 5	9,80	5			10

15 25 9	3,05	00	25	3,0	10
15 25 0	3,05	0			10
15 25 1	3,05	1			10
15 25 2	3,05	2			10
15 25 3	3,05	3			10
15 25 4	3,50	4			10
15 25 5	13,30	5			10


Best.-Nr.	€/Stück	Por.	B	*S	VE
Cat. -No.	€/piece		mm	ca. mm	
15 30 9	3,35	00	30	3,8	10
15 30 0	3,35	0			10
15 30 1	3,35	1			10
15 30 2	3,35	2			10
15 30 3	3,35	3			10
15 30 4	4,05	4			10
15 30 5	14,40	5			10
15 40 9	4,45	00	40	4,5	10
15 40 0	4,45	0			10
15 40 1	4,45	1			10
15 40 2	4,45	2			10
15 40 3	4,45	3			10
15 40 4	5,30	4			10
15 40 5	18,65	5			10
15 50 9	5,75	00	50	4,5	10
15 50 0	5,75	0			10
15 50 1	5,75	1			10
15 50 2	5,75	2			10
15 50 3	5,75	3			10
15 50 4	6,75	4			10
15 50 5	20,95	5			10
15 60 9	9,35	00	60	5,0	5
15 60 0	9,35	0			5
15 60 1	9,35	1			5
15 60 2	9,35	2			5
15 60 3	9,35	3			5
15 60 4	11,75	4			5
15 60 5	37,45	5			5
15 65 9	13,60	00	65	5,5	5
15 65 0	13,60	0			5
15 65 1	13,60	1			5
15 65 2	13,60	2			5
15 65 3	13,60	3			5
15 65 4	16,45	4			5
15 65 5	39,50	5			5
15 70 9	16,25	00	70	5,5	5
15 70 0	16,25	0			5
15 70 1	16,25	1			5
15 70 2	16,25	2			5
15 70 3	16,25	3			5
15 70 4	21,10	4			5
15 70 5	47,20	5			5
15 80 9	17,65	00	80	6,2	5
15 80 0	17,65	0			5
15 80 1	17,65	1			5
15 80 2	17,65	2			5
15 80 3	17,65	3			5
15 80 4	22,90	4			5
15 80 5	57,70	5			5
15 90 9	20,95	00	90	7,0	5
15 90 0	20,95	0			5
15 90 1	20,95	1			5
15 90 2	20,95	2			5
15 90 3	20,95	3			5
15 90 4	27,50	4			5
15 90 5	64,15	5			5

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
**Please contact us for your special designs and intermediate sizes**

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
Please contact us for your special designs and intermediate sizes

Preise in €uro netto/Stück  
Prices in €uro currency net/piece

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
**Filter discs**

edge not fused, centred, bi-plane ground  
The faces of our filterdiscs from the 16 and 61 – 63 series are precisely ground down to tight tolerances and plane surfaces as may be required in specific installation circumstances.

Best. -Nr. Cat. -No.	€/Stück €/piece	Por.	B mm	S mm	VE
16 60 9	16,05	00	60	4,0	5
16 60 0	16,05	0			5
16 60 1	16,05	1			5
16 60 2	16,05	2			5
16 60 3	16,05	3			5
16 60 4	19,20	4			5
16 60 5	44,90	5			5
16 70 9	24,75	00	70	4,5	5
16 70 0	24,75	0			5
16 70 1	24,75	1			5
16 70 2	24,75	2			5
16 70 3	24,75	3			5
16 70 4	29,25	4			5
16 70 5	55,35	5			5
16 80 9	29,05	00	80	5,0	5
16 80 0	29,05	0			5
16 80 1	29,05	1			5
16 80 2	29,05	2			5
16 80 3	29,05	3			5
16 80 4	34,05	4			5
16 80 5	68,85	5			5
16 90 9	36,95	00	90	5,5	5
16 90 0	36,95	0			5
16 90 1	36,95	1			5
16 90 2	36,95	2			5
16 90 3	36,95	3			5
16 90 4	43,00	4			5
16 90 5	79,65	5			5
16 11 9	61,75	00	100	6,8	1
16 11 0	61,75	0			1
16 11 1	61,75	1			1
16 11 2	61,75	2			1
16 11 3	61,75	3			1
16 11 4	70,50	4			1
16 11 5	128,50	5			1
16 12 9	85,50	00	120	7,5	1
16 12 0	85,50	0			1
16 12 1	85,50	1			1
16 12 2	85,50	2			1
16 12 3	85,50	3			1
16 12 4	97,50	4			1
16 12 5	158,85	5			1

Best. -Nr. Cat. -No.	€/Stück €/piece	Por.	B mm	S mm	VE
16 10 9	4,50	00	10	2,0	10
16 10 0	4,50	0			10
16 10 1	4,50	1			10
16 10 2	4,50	2			10
16 10 3	4,50	3			10
16 10 4	5,05	4			10
16 10 5	11,40	5			10
16 20 9	4,80	00	20	2,5	10
16 20 0	4,80	0			10
16 20 1	4,80	1			10
16 20 2	4,80	2			10
16 20 3	4,80	3			10
16 20 4	5,40	4			10
16 20 5	11,85	5			10
16 25 9	5,05	00	25	2,5	10
16 25 0	5,05	0			10
16 25 1	5,05	1			10
16 25 2	5,05	2			10
16 25 3	5,05	3			10
16 25 4	5,60	4			10
16 25 5	15,40	5			10
16 30 9	5,35	00	30	3,0	10
16 30 0	5,35	0			10
16 30 1	5,35	1			10
16 30 2	5,35	2			10
16 30 3	5,35	3			10
16 30 4	5,90	4			10
16 30 5	16,25	5			10
16 40 9	7,20	00	40	3,5	10
16 40 0	7,20	0			10
16 40 1	7,20	1			10
16 40 2	7,20	2			10
16 40 3	7,20	3			10
16 40 4	8,00	4			10
16 40 5	21,35	5			10
16 50 9	8,80	00	50	3,8	10
16 50 0	8,80	0			10
16 50 1	8,80	1			10
16 50 2	8,80	2			10
16 50 3	8,80	3			10
16 50 4	9,90	4			10
16 50 5	24,10	5			10

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Preise in Euro netto/Stück  
Prices in Euro currency net/piece

**Please contact us for any special designs.**

Reproduction of this catalogue, in whole or in part, is permitted only with our written consent.

Technical information, dimensions and drawings in this catalogue are only of approximate value. Deviations are possible and all items are subject to technical changes.

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**Introduction**

Borosilicate glass 3.3 is specified as a standard glass - it contains mainly sand, calcium carbonate and sodium carbonate.

High quality sand, dehydrated borate boric acid, alumina and salt are used in its manufacture and strict control is applied to ensure the specifications for purity and consistency of these materials are maintained. Because of its resistance to heat, corrosion and thermal shock borosilicate glass 3.3 is used extensively in science and in industry.

**Fusion of Sintered Glass**

All **VitraPOR®** sintered filters in our glass apparatus are solely joint to the massive glass by fusing, i.e. welding without additives. No foreign substances like adhesives impair the inert properties of the pure borosilicate glass 3.3.

**Standard Thicknesses Filterdiscs (Art. 15 - ...)**

Please note, that the values are only of approximate value and given without tolerances

Outer Diameter	Approximate Thickness
Ø 5 – 10 mm	2,8 mm
Ø 20 – 25 mm	3,0 mm
Ø 30 mm	3,8 mm
Ø 40 – 50 mm	4,5 mm
Ø 60 mm	5,0 mm
Ø 65 – 70 mm	5,5 mm
Ø 80 mm	6,2 mm
Ø 90 mm	7,0 mm
Ø 100 – 120 mm	8,5 mm
Ø 130 – 150 mm	10,0 mm
Ø 160 – 180 mm	13,0 mm
Ø 190 – 200 mm	15,0 mm
Ø 210 – 250 mm	19,0 mm
Ø 260 – 300 mm	26,0 mm
Ø 350 mm	28,0 mm
Ø 380 mm	38,0 mm

The material-specific shrinkage of finer pore sizes results in correspondingly smaller thicknesses of the filters.

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at 20°C 1.13 W/m·K  
150°C 1.25 W/m·K  
300°C 1.42 W/m·K

**Chemical Composition**

The following is a typical analysis of borosilicate glass 3.3:

Element	% by weight
Silica (SiO <sub>2</sub> )	80.60
Boric oxide (B <sub>2</sub> O <sub>3</sub> )	12.60
Sodium oxide (Na <sub>2</sub> O)	4.20
Alumina (Al <sub>2</sub> O <sub>3</sub> )	2.20
Iron oxide (Fe <sub>2</sub> O <sub>3</sub> )	0.04
Calcium oxide (CaO)	0.10
Magnesium oxide (MgO)	0.05
Chlorine (Cl)	0.10

There may be heavy metals in a concentration of less than 5 ppm.

**Physical Properties**

Coefficient of Expansion	33 x 10 <sup>-7</sup> /K	between 20°C - 300°C
Specific Heat	0,8x10 <sup>3</sup> J/kg K <sup>-1</sup>	at 20°C
Thermal Conductivity	1.13	W/m K <sup>-1</sup> at 20°C
Densit	2.23 x	10 <sup>3</sup> kg/m <sup>3</sup>
Poisson's Ratio	0.22	between 25°C - 400°C
Young's Modulus	6.500	kg/mm <sup>2</sup> at 25°C
Rigidity Modulus E	64 x 10 <sup>3</sup>	MPa
DPH (Vickers) Hardness	580	Kg/mm <sup>2</sup> with 50 gram load
Relative Hardness	1.52	(comparative Soda-Lime = 1.0)
Refractive Index	1.474	Sodium D - Line
Dielectric Constant	4.6	at 1 MHz and 20°C
Loss Factor	2.6	% at 1 MHz and 20°C
Log <sub>10</sub> Volume Resistivity	15	Ohm · cm/s at 20°C
Surface Resistivity	10 <sup>13</sup>	Ohm s·cm <sup>-2</sup> at 50% humidity
Transformation Point (t <sub>g</sub> )	525	°C

**Thermal Expansion**

The coefficient of linear thermal expansion is one of the characteristic properties by which glass is measured and is defined as the change in unit length per degree rise in temperature. The thermal expansion is 33 x 10<sup>-7</sup> per °C between 20°C and 300°C and is controlled to a tolerance of ± 0.5 x 10<sup>-7</sup> per °C. Due to its porous structure, sintered glass has a slightly different thermal expansion. Care should be taken when fusing it to massive borosilicate glass 3.3.

**Thermal Resistance**

Since borosilicate glass has a high softening point, it can be used at temperatures of the order of 515°C. Care must be exercised, however, in cooling down apparatus which has been held at high temperatures. The cooling process must be slow and uniform, particularly in the case of thick sections (see annealing). If a temperature gradient is applied to the glass, stresses are set up - the hotter glass being under compression and the cooler under tension. The magnitude of the tensile stress depends on the temperature difference which in turn depends on the glass thickness.

ATTENTION:

Abrasions greatly reduce the temperature shock resistance.

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
## Maximum Temperatures

In general it is recommended, that the strain point (515°C) be regarded as the maximum safe operating temperature of borosilicate glassware 3.3. For some shapes and for a relatively short time this limit can be exceeded, but at 580°C there is danger of deformation and in case of sintered glassware the porous structure may be changed. At high temperatures the glass may acquire permanent stress on cooling and this may result in subsequent breakage.

If it is suspected that permanent stress has occurred, the article should be annealed making references to the annealing process explained below. Permanent stress can greatly reduce the mechanical and thermal resistance.

## Viscosity

Although glass has no defined softening point, there are four points in the viscosity-/temperature relation which have accepted definitions (ISO 7884-2/-3/-4):

1. Strain Point	515°C		
Viscosity of the glass is		$10^{14.5}$	dPa·s
2. Annealing Point	565°C		
Viscosity of the glass is		$10^{13.0}$	dPa·s
3. Softening Point	820°C		
Viscosity of the glass is		$10^{7.6}$	dPa·s
4. Working Point	1250°C		
Viscosity of the glass is		$10^{4.0}$	dPa·s

## Annealing

The annealing of glass is the process by which it is heated to and held at a controlled temperature for a defined period to relieve internal stresses. Careful cooling under controlled conditions is essential to ensure, that no stresses are reintroduced by chilling.

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### Tabelle / Table

Refractive index (sodium D line)	1.474
Dispersion (NF-NC)	0.00738
Visible light transmission	92% in 2 mm thick glass 91% in 6 mm thick glass
Stress optical coefficient	3.5 nm/kg cm <sup>2</sup> per cm length

## Chemical Resistance

The chemical resistance of glass is better than that of all other known materials. It shows resistance to water, acids and alkalis, salt, organic substances as well as chlorine and bromine. Hydrofluoric acid, concentrated phosphoric acid and strong alkalis attack the glass surface at higher concentrations and temperatures.

The following are typical results of tests undertaken to international standards:

### 1. Hydrolytic resistance

Glass grains (300 - 500 microns)  
acc. ISO 719 (98°C) Water attack: < 0,01 mg Na<sub>2</sub>O/g grains  
Class ISO 719 - HGB 1  
acc. ISO 720 (121°C)  
Class ISO 720 - HGA 1

### 2. Acid Resistance

Resistance acc. DIN 12116 Acid Class 1  
Resistance acc. DIN/ISO 1766 = Loss in weight < 100 µg Na<sub>2</sub>O/dm<sup>2</sup>

### 3. Alkali Resistance

Resistance acc. DIN 52322 / ISO 695 Class A2

## Optical Properties

Massive borosilicate glass 3.3 has a high ability to transmit light throughout the visible range of the spectrum. Ranging from approx. 310 to 2300 nm there is virtually no absorption through the thickness of glass generally encountered in laboratory ware (up to approx. 6 mm). The following table shows additional data:

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






**Porosity**

The VitraPOR glass filters are divided into seven porosity classes from 1,0 µm to 500 µm. Therefore they are suitable for all applications in the chemical and pharmaceutical industry, in laboratories as well as in large technical plants.

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
Druckunterschied/Pressure differential  $\Delta p$  (mbar)

## Flow rate and Pressure Drop


Flow rates of liquids and drop in pressure of gases are functions of:

- diameter of the filter
- porosity of the filter
- relation of pressure/vacuum in the apparatus
- physical properties of the medium to be filtered.

To determine the possible application of glass filter-products and to choose the best apparatus for the work intended, it is necessary to know the flow rates of liquids and gases. Please use our diagrams 1 and 2 which are shown for water and air. The data apply to filter discs of 30 mm diameter with average properties and usually effective filter area and - thicknesses, as well as poresizes.

The flow rates for other diameters can easily be taken by multiplying this factor by the conversion factor. Minor deviations with very limited tolerances are possible. Strict control is recommended, especially when filtering bacteria using porosity 5.

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## Cleaning

Successful experimental results can only be achieved by using clean apparatus. The following advice may be used for the massive and porous borosilicate glass 3.3:

Wash glassware as quickly as possible after use. During washing with water all parts of the article should be scrubbed with a soft brush. Avoid any abrasion of the glassware! After cleaning thoroughly rinse with distilled water. Drying can be undertaken either in air or at a temperature not exceeding 100°C.

If after mechanical cleaning for example filter-pores should still remain clogged, the thorough chemical cleaning is required. A few suggestions found generally useful follow:

Contamination	Solvent
<b>Barium sulphate</b>	Hot concentrated sulphuric acid (100°C), Silver-chloride Ammonian or Sodium hyposulfite
<b>Copper or Iron Oxides</b>	Hot Hydrochloric acid plus potassium chlorate
<b>Mercury residue</b>	Hot Nitric acid
<b>Mercury sulphide</b>	Hot aqua regia
<b>Albumen</b>	Hot hydrochloric acid or hot ammonia
<b>Grease, oil, fatty materials</b>	Carbon tetrachloride
<b>Organic matters</b>	Hot concentrated cleaning solution, or hot concentrated sulfuric acid plus a few drops of sodium or potassium nitrite
<b>Glucose</b>	Carefully heat to approx. 200°C with mixed acid (5% H <sub>2</sub> SO <sub>4</sub> and 1% HNO <sub>3</sub> )

Extensive rinsing with water must be obviously follow.

### Warning!

Hot concentrated phosphoric acid and hot alkaline solutions attack the glass surface, they are unsuitable as cleaning agents. If they have to be filtered, an increase in pore size and reduced life of the fritted disc is unavoidable.

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## Safety Instructions

### Cleaning

Before the first use of glass-filters, they must be cleaned with diluted, hot hydrochloric acid, followed by several rinses with distilled water. Glass-filters should be cleaned immediately after their use. The surfaces can be wiped with a squeegee or a brush. Surface rinsing with a wash bottle or back-flushing under the water tap is often sufficient.

### Pressure and Vacuum

Due to the mechanical properties, a vacuum- or pressure resistance of glass cannot be guaranteed. Scratches or internal abrasions of the glass surface can greatly reduce the mechanical and thermal strength. Filter pores may clog during filtration and cause a pressure drop - in any case 0,1 MPa may not be exceeded. Never subject glassware to sudden pressure changes.

### Temperature Resistance

The maximum allowable service temperature is approx. 500°C. To avoid permanent stress, careful heating and cooling according to the schedule for borosilicate glass 3.3 should be ensured. Permanent stress can reduce the mechanical and thermal resistance. Never subject glassware to sudden temperature changes.

### Chemical resistance

High concentrations of hydrofluoric acid or hot alkali solutions attack the glass and can cause deterioration in the filter and an increase of the pore size.

### Safety advice

When working with glassware always wear protective glasses and protective gloves to avoid injuries. Use a safety screen, -hood or similar protection when working with glassware subject to pressure or vacuum. Be prepared, that reagents may leak from a broken vessel.

Please heed these guidelines in combination with the respective state-specific regulations for the use of glassware in the laboratory.

All information also available for download:  
[www.robuglas.com](http://www.robuglas.com)

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

## Twin-Filter

The twin-filter combines two different porosities sintered to a single filter. A coarser layer serves as a pre-filter and prevents the finer pores from clogging. The finer layer holds back even smaller particles and guarantees an optimal filtration result. This multi-layer technology extends the life of the filter and increases flowrates and mechanical stability.

## Casting Technologies

Advanced casting technologies allow the moulding of sintered glass to a large variety of shapes and sizes with porous or even massive structure. Pore diameters of below 1,0 micron at porosities of 35% can be achieved.

## Porous Glass Beads

ROBU® has developed a new carrier substrate for cell cultivation from the proven family of VitraPOR® sintered filter products. Pure, binder-free USP-class I borosilicate glass 3.3 is manufactured by our proprietary process to form 1.5 mm or alternatively 4 mm diameter beads with defined pore sizes of approximately 60 micrometers. The special process uses no corrosives and produces a porous structure with smooth, yet large surface areas. The glass is completely bio-compatible and its high electrical surface charge (Zeta-Potential) helps to stimulate cell growth. As a result, you can yield more than 70% viable cells from your cultivation process - far superior to conventional PMMA and PVC carriers.

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## Terms of delivery and payment

1. The following terms of delivery and payment apply to all consignments made and all services supplied by us. In case purchasing conditions of our customers depart from our terms, they shall apply only after having been expressly confirmed by us in writing. The acceptance of goods supplied shall be deemed as recognition of our terms.

2. Prices charged for are those in effect on the day of delivery and, unless otherwise provided, are to be understood in EURO currency ex works plus packing, dispatch and insurance at the expense of our customers, payable immediately net cash without discount. At the time of shipment, risk passes to the buyer as soon as we have handed over the goods to the forwarding agents chosen by us.

2.1 The return of packaging material will not be accepted, unless we are obliged to do so by German laws and regulations.

3. In case payment dates are exceeded, all legal consequences of default may take effect without any prior formal notice. We reserve the right to charge interest at the rate applied from time to time by our bankers. The total balance will be due for immediate payment regardless of any other terms of payment provided. If there will be doubt in creditability, even if the order is already reconfirmed, we have the right to ask for prepayment of irrevocable letter of credit or to terminate the contract.

4. For special designs we reserve the right to lift prices within reasonable limits as well as to adequately vary the amounts agreed upon for delivery; surplus quantities have to be accepted. Related to quantities ordered a +/- 10% deviation is allowed. Special designs cannot be returned and orders not be cancelled without our written consent. Moulds for special designs ordered are due to the account of the customer, but still remain our property.

5. We do our best to keep the delivery terms stated. Due to the risk and peculiarities of glassworking, these terms are non committal, unless otherwise expressly stipulated. Our contractual obligations are subject to our being supplied by our subcontractors correctly and in due course of time.

6. Should there be any reason for complaints, these shall be brought to formal notice within eight days after receipt of goods, otherwise the goods shall be deemed as approved. We may choose either to repair the defect or to supply substitute goods. In case of any damages, including consequential damages, we shall be liable at the most up to a value not exceeding the amount invoiced. Parts subject to wear and tear are excluded from warranty for defects. To proof any complaints, the delivered goods have to be returned.

6.1 Claims on the part of infringement of proprietary rights are excluded, if the infringement is a result of special instructions or a not foreseeable application or use of the goods by the Customer.

7. Drawings, photographs, dimensions and weights are only of approximate value, unless they have expressly been stated as mandatory. The Buyer shall be responsible for working drawings submitted by him not to encroach upon patent rights of third parties. He shall save us harmless against claims for compensation. Information submitted shall not be made accessible to third parties or to be reproduced or to be used for purposes other than agreed upon.

8. Goods supplied will remain our property until all of our claims, including also those arising in the future, have been paid in full. The Buyer is entitled to process the goods received and to re-sell them, provided the following conditions will be complied with:

8.1 In case goods are processed or transformed by the Buyer, we are deemed, within the meaning of Article 950 of the German Code of Civil Law (BGB), the manufacturer and acquire title of ownership to intermediate or finished products. The processing party is only a depositary. If goods reserved are incorporated into or processed together with other items not owned by us, we acquire title of co-ownership to the new object in proportion of the value of the goods reserved to the value of said other items.

8.2 Goods are allowed to be sold only by way of ordinary and regular business transactions and only if claims resulting from re-sales have not previously been assigned to third parties. Claims resulting from re-sales, and to which the Buyer is entitled, shall be considered as assigned to us from the moment the buyer enters into a sales agreement with us, and also in as far as our goods are incorporated into or processed together with other items. In said case, the claims assigned serve as a guarantee for us only to an amount not exceeding the sales price of goods reserved. We will not collect the claims assigned, as long as the Buyer meets his liabilities. The Buyer shall be obliged to name to us, on request, the garnishes and to notify to these said assignments. He shall be entitled to collect claims himself, as long as he shall not receive from instructions to the contrary. Amounts collected by him shall immediately be paid to us to the extent our claims are due.

8.3 Mortgages or chattel mortgages of goods reserved or of claims assigned are not allowed. The Buyer shall notify us immediately of eventual seizures by third parties of goods supplied with retention of title or of claims assigned.

9. Should, in event of sales to foreign countries, the retention of title not be admitted with the same effect as under German law, the goods involved will remain our property until payment of all our claims resulting from the contractual relationship brought about by the sale of said goods. In case this retention of title is neither admitted with the same effect as under German law but if it is allowed to reserve other rights to said goods, we will be entitled to exercise all these rights. The Buyer shall be obliged to cooperate in our measures we intend to take for the protection of our right of ownership or, in its place, of any other title of the goods.

10. All sale contracts are due to the law of the Federal Republic of Germany. Any law case as well as matters concerning bill of exchange, drafts, etc. have to be placed with the competent district court of our company. In case we appear plaintiff, we will be entitled to enter an action also at the purchasers place of business.

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With the background of a long term experience **ROBU®** has developed a new and flexible porous quartz membrane. It is produced in a special casting process using a specialized technology. It not only offers various advantages in electrochemical, opto-electronic and electrokinetic applications, but also in the medical field or as a filter and support structure:

- Raw material pure, amorphous SiO<sub>2</sub>
- Temperature resistant to more than 1.000 °C
- Superior temperature change resistance
- Very good chemical resistance
- High electric potential, porosity and large surface area
- Mono- or multimodal pore sizes from 0.5 to over 5 µm dia.
- Tape thicknesses from below 1.0 to over 10 mm

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Our wide range of porosities and sizes makes it easy for you. If your specific demands cannot be met with our standard products, we will work closely with you to create what you need. ROBU can fabricate sinterfilters in almost any shape and size.

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Our deliveries are fast and reliable! Most standard products and many specialty items are available within a few days.

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Our DIN/ISO 9001 certified quality management and in process controls of all major properties of the sintered elements ensure the consistent quality of all products.

**ROBU**<sup>®</sup>

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Schützenstr. 13, D-57644 Hattert, Germany

Tel.: ..49-(0) 2662/8004-0

Fax.: ..49-(0) 2662/8004-40

[www.robuglas.com](http://www.robuglas.com)

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ROBU Glasfilter-Geräte GmbH