


Urine Analysis

Combi 10[®] SGL Optimized Urine Test Strip



Best available vitamin C protection

+ Safe, early diagnosis of serious diseases

Desiccant enclosed in cap

+ Convenient and safe

Premium color charts perfectly match reaction colors



+ Clear read-off, reliable results

MACHEREY-NAGEL

www.mn-net.com




Since 1911



Best available vitamin C protection

Comparative data clearly show: Combi 10® SGL urine test strips have the best protection of blood and glucose tests against vitamin C. With other brands, vitamin C may prevent the oxidation reaction essential for these test. This may cause false negative results for blood and glucose. Using Medi-Test Combi 10® SGL, **blood and glucose** in urine **will be safely detected**. Patients can receive treatment at the earliest stage thus increasing chances for success.




► Concentration of vitamin C in the sample with no effect on the results for blood and glucose in urine (mg/dL)



Desiccant enclosed in cap

Test strips Medi-Test Combi 10® SGL come in tubes with a desiccant containing cap. The desiccant ensures a low humidity in the tube and the longevity of the test strips. Many other packages contain a simple desiccant pack that may easily be spoiled or accidentally removed. The desiccant in the cap is **convenient and safe**. Closing the cap safely protects the test strips.



High quality color chart

The reaction colors on the test strip exactly match the colors on the label. This is ensured by a high quality printing process using specially mixed full colors. Other color charts often yield far less precise matches that may lead to inaccurate results. The evaluation of Medi-Test Combi 10® SGL is always **clear and reliable** providing perfectly accurate results.

Your local distributor:



KATEN100032/FlyerCombi10SGL_en4/13/2.2013 DD
Printed in Germany

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

Germany and international:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

Switzerland:
MACHEREY-NAGEL AG
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

France:
MACHEREY-NAGEL EURL
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

USA:
MACHEREY-NAGEL Inc.
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com

