


## Couloir à manche inox 304



Panneau repoussé . Fabrication mécano soudée

| Réf : | A   | a   | D   | Volume |
|-------|-----|-----|-----|--------|
| 8800  | 190 | 100 | 50  | 0,12 L |
| 8801  | 250 | 145 | 75  | 0,45 L |
| 8802  | 320 | 190 | 105 | 1,00 L |
| 8803  | 350 | 230 | 120 | 1,50 L |
| 8804  | 370 | 250 | 150 | 2,00 L |

## Pelle inox 304



| Réf :    | A   | D   | m    | Poids | Matière  |
|----------|-----|-----|------|-------|----------|
| 343-250  | 180 | 60  | 250  | 80g   | Inox 304 |
| 343-500  | 205 | 75  | 500  | 110g  | Inox 304 |
| 343-750  | 230 | 90  | 750  | 140g  | Inox 304 |
| 343-1000 | 270 | 110 | 1000 | 200g  | Inox 304 |

## Pelles, aluminium



| Réf : | A       | a       | D  | ø | m    |
|-------|---------|---------|----|---|------|
|       | longeur | longeur |    |   |      |
| 8814  | 65      | 40      | 12 |   | 5    |
| 8815  | 85      | 52      | 13 |   | 10   |
| 8816  | 135     | 80      | 24 |   | 50   |
| 8817  | 155     | 95      | 30 |   | 65   |
| 8818  | 185     | 115     | 32 |   | 90   |
| 8819  | 210     | 120     | 37 |   | 160  |
| 8820  | 250     | 145     | 45 |   | 240  |
| 8821  | 290     | 175     | 55 |   | 400  |
| 8822  | 310     | 185     | 57 |   | 480  |
| 8823  | 350     | 220     | 63 |   | 900  |
| 8824  | 400     | 270     | 83 |   | 1600 |
| 8825  | 430     | 300     | 95 |   | 2500 |
| 8827  | 250     | 120     | 35 |   | 90   |
| 8828  | 285     | 120     | 35 |   | 90   |
| 8829  | 365     | 120     | 35 |   | 90   |

## Pelles pour balance



| Réf.  | A        | D/ø    | Matière   |
|-------|----------|--------|-----------|
|       | Longueur |        |           |
| 12800 | 70       | 32     | INOX 304  |
| 12801 | 90       | 32     | INOX 304  |
| 12802 | 105      | 32     | INOX 304  |
| 12803 | 120      | 32     | INOX 304  |
| 12804 | 120      | 30-65  | INOX 304  |
| 12805 | 160      | 30-75  | INOX 304  |
| 12806 | 200      | 55-100 | INOX 304  |
| 12810 | 70       | 32     | Aluminium |
| 12811 | 90       | 32     | Aluminium |
| 12812 | 105      | 32     | Aluminium |
| 12813 | 120      | 32     | Aluminium |

| Réf.     | prix HT € | Réf.  | prix HT € |
|----------|-----------|-------|-----------|
| 8800     | 24,67     | 8814  | 1,88      |
| 8801     | 27,18     | 8815  | 2,22      |
| 8802     | 32,12     | 8816  | 4,31      |
| 8803     | 37,02     | 8817  | 5,44      |
| 8804     | 51,87     | 8818  | 7,42      |
| 343-250  | 16,83     | 8819  | 9,42      |
| 343-500  | 21,64     | 8820  | 11,98     |
| 343-750  | 27,24     | 8821  | 18,31     |
| 343-1000 | 32,56     | 8822  | 21,53     |
| 12800    | 16,11     | 8823  | 31,40     |
| 12801    | 16,69     | 8824  | 52,85     |
| 12802    | 17,25     | 8825  | 71,70     |
| 12803    | 17,93     | 8827  | 14,20     |
| 12804    | 52,76     | 8828  | 15,86     |
| 12805    | 57,92     | 8829  | 16,86     |
| 12806    | 68,76     | 12812 | 10,73     |
| 12810    | 10,16     | 12813 | 11,64     |
| 12811    | 10,06     |       |           |