





FABRICATION SPECIALE
SUR DEMANDE

DIMENSIONNEMENT




CHOIX DE LA MATIERE DU TUYAU :

- **Silicone USP6 spéciale pharma**
Renfort textile ou
Spirale noyée dans le silicone



- **PTFE renforcé tresse métallique, spécial chimie**





*FABRICATION SPECIALE
SUR DEMANDE*

CHOIX DU RACCORD :


- Par embouts sertis en inox 316L monoblocs sans soudure suivant le système **NRS®** (Non Retention System)
- Avec le raccord inox breveté **MSP®** (Montage sur place)
- Tout type de raccordement possible : SMS, DIN, CLAMP, RJT, ISS, FIL, GAZ, NPT, JIC, BRIDES, etc...



Raccord MSP® breveté
(Montage sur place)
*MSP® patented fittings
(Assembly on site)*
Disponible du Ø 13 au Ø 51
Available from Ø13 to Ø 51

® eingetragenes Markenzeichen GECI - ® GECI registered trade mark

- *The 316L stainless steel swaged fittings are monoblocs (without welding), in accordance with the NRS® design (Non Retention System).*
- *Stainless steel fitting with patented MSP® design. Fittings are reusable.*
- *Many kinds of fittings: SMS, DIN, CLAMP, RJT, ISS, FIL, GAZ, NPT, JIC, FLANGES, etc...*



Concept NRS®
(Non Retention System)
*NRS® Design
(Non Retention System)*

MARQUAGE – IDENTIFICATION – QUALITE :

Flexibles et raccords identifiés unitairement

Bracelet de couleur en option

certificats FDA, USP6, et électro polissage sur demande.