

Magnetic Stirrers | Accessories



Synthesis blocks

Dry heating block square series

	Ident. No.	Price
H 135.10 Square carrier with handle	0004448200	on request
H 135.10 Square carrier without handle	0025000832	94,- EUR
H 135.101 Block 16 x 4 ml	0025000626	35,- EUR
H 135.102 Block 16 x 8 ml	0025000627	38,50 EUR
H 135.103 Block 9 x 16 ml	0025000628	38,50 EUR
H 135.104 Block 4 x 20 ml	0025000629	33,- EUR
H 135.105 Block 4 x 30 ml	0025000630	38,50 EUR
H 135.106 Block 4 x 40 ml	0025000631	42,50 EUR
H 135.107 Block 100 ml	0025000632	34,- EUR
H 135.108 Block 250 ml	0025000633	40,50 EUR



100 ml Flask heating block series

	Ident. No.	Price
H 135.20 Flask carrier 100 ml with handle	0004448300	on request
H 135.20 Flask carrier 100 ml without handle	0025000634	82,50 EUR
H 135.201 Flask inlay 10 ml	0025000636	19,80 EUR
H 135.202 Flask inlay 25 ml	0025000637	19,80 EUR
H 135.203 Flask inlay 50 ml	0025000638	19,80 EUR



500 ml Flask heating block series

	Ident. No.	Price
H 135.30 Flask carrier 500 ml with handle	0004448400	on request
H 135.30 Flask carrier 500 ml without handle	0025000639	94,- EUR
H 135.301 Flask inlay 100 ml	0025000641	44,- EUR
H 135.302 Flask inlay 250 ml	0025000642	45,- EUR



1000 ml Flask heating block series

	Ident. No.	Price
H 135.40 Flask carrier 1000 ml with handle	0004448500	on request
H 135.40 Flask carrier 1000 ml without handle	0025000833	136,- EUR
H 135.401 Flask inlay 500 ml	0025000644	77,- EUR



2000 ml Flask heating block series

	Ident. No.	Price
H 135.50 Flask carrier 2000 ml with handle	0004448600	on request
H 135.50 Flask carrier 2000 ml without handle	0025000834	168,- EUR
H 135.501 Flask inlay 1000 ml	0025000645	113,- EUR



3000 ml Flask heating block series

	Ident. No.	Price
H 135.60 Flask carrier 3000 ml with handle	0004448700	on request
H 135.60 Flask carrier 3000 ml without handle	0025000835	244,- EUR
H 135.601 Flask inlay 2000 ml	0025000648	179,- EUR



5000 ml Flask heating block series

	Ident. No.	Price
H 135.70 Flask carrier 5000 ml with handle	0004448800	on request
H 135.70 Flask carrier 5000 ml without handle	0025000836	367,- EUR
H 135.701 Flask inlay 3000 ml	0025000651	313,- EUR

AVAILABLE 2014
(Items in green colored text)



www.ika.com



www.facebook.com/
IKAworlwide



www.twitter.com/
IKAworlwide



www.instagram.com/
IKAworlwide

Stirring bars



IKAFLON® round, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 10 Set (5 Pcs)	10 x 6 mm	0004488600	6,10 EUR
IKAFLON® 15 Set (5 Pcs)	15 x 6 mm	0004488700	6,10 EUR
IKAFLON® 20 Set (5 Pcs)	20 x 8 mm	0004488800	6,10 EUR
IKAFLON® 25 Set (5 Pcs)	25 x 8 mm	0004488900	9,20 EUR
IKAFLON® 30 Set (5 Pcs)	30 x 8 mm	0004489000	7,70 EUR
IKAFLON® 40 Set (5 Pcs)	40 x 8 mm	0004489100	9,30 EUR
IKAFLON® 50 Set (5 Pcs)	50 x 8 mm	0004489200	11,10 EUR
IKAFLON® 80 Set (5 Pcs)	80 x 10 mm	0004489300	26,50 EUR




IKAFLON® glass round

	(L x Ø)	Ident. No.	Price
IKAFLON® 25 Set (5 Pcs)	25 x 8 mm	0004492200	10,80 EUR
IKAFLON® 30 Set (5 Pcs)	30 x 8 mm	0004492400	11,30 EUR
IKAFLON® 40 Set (5 Pcs)	40 x 8 mm	0004492600	12,30 EUR
IKAFLON® 50 Set (5 Pcs)	50 x 8 mm	0004492800	13,30 EUR



IKAFLON® power SmSo, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 15 Set (5 Pcs)	15 x 9 mm	0004493000	17,60 EUR
IKAFLON® 30 Set (5 Pcs)	30 x 12 mm	0004493200	26,50 EUR
IKAFLON® 50 Set (5 Pcs)	50 x 21 mm	0004493400	32,50 EUR




IKAFLON® slide round, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 25 Set (5 Pcs)	25 x 6 mm	0004493800	7,60 EUR
IKAFLON® 30 Set (5 Pcs)	30 x 6 mm	0004494000	8,20 EUR
IKAFLON® 40 Set (5 Pcs)	40 x 8 mm	0004494200	13,60 EUR
IKAFLON® 50 Set (5 Pcs)	50 x 8 mm	0004494400	15,70 EUR



IKAFLON® crown, PTFE-coated

	(Ø x H, for beaker)	Ident. No.	Price
IKAFLON® 9 Set (5 Pcs)	9 x 6 mm	0004498400	8,50 EUR



IKAFLON® ellipse, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 20 Set (5 Pcs)	20 x 10 mm	0004494600	12,- EUR
IKAFLON® 25 Set (5 Pcs)	25 x 12 mm	0004494800	13,20 EUR
IKAFLON® 30 Set (5 Pcs)	30 x 15 mm	0004495000	17,30 EUR
IKAFLON® 35 Set (5 Pcs)	35 x 15 mm	0004495200	22,50 EUR
IKAFLON® 40 Set (5 Pcs)	40 x 20 mm	0004495400	28,- EUR
IKAFLON® 50 Set (5 Pcs)	50 x 20 mm	0004495600	39,50 EUR
IKAFLON® 60 Set (5 Pcs)	60 x 20 mm	0004495800	57,50 EUR
IKAFLON® 70 Set (5 Pcs)	70 x 20 mm	0004496000	69,- EUR



IKAFLON® cross, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 10 Set (5 Pcs)	10 x 10 mm	0004496200	10,30 EUR
IKAFLON® 20 Set (5 Pcs)	20 x 20 mm	0004496400	13,20 EUR
IKAFLON® 25 Set (5 Pcs)	25 x 25 mm	0004496600	15,- EUR
IKAFLON® 30 Set (5 Pcs)	30 x 30 mm	0004496800	19,50 EUR
IKAFLON® 38 Set (5 Pcs)	38 x 38 mm	0004497000	24,50 EUR




IKAFLON® bone, PTFE-coated

	(L x Ø)	Ident. No.	Price
IKAFLON® 37 Set (5 Pcs)	37 x 8 x 20 mm	0004497200	35,50 EUR
IKAFLON® 54 Set (5 Pcs)	54 x 8 x 20 mm	0004497400	42,50 EUR



TRIKA® PTFE-coated

	(L x Ø)	Ident. No.	Price
TRIKA® 25 Set (5 Pcs)	25 x 8 mm	0004499300	8,30 EUR
TRIKA® 40 Set (5 Pcs)	40 x 8 mm	0004499400	12,80 EUR
TRIKA® 55 Set (5 Pcs)	55 x 14 mm	0004499500	20,50 EUR
TRIKA® 80 Set (5 Pcs)	80 x 18 mm	0004499600	28,50 EUR



IKAFLON® beaker, PTFE-coated

	(Ø x H, for beaker)	Ident. No.	Price
IKAFLON® 67	67 x 21 mm with 50 x 8	0004497600	197,- EUR
IKAFLON® 74	74 x 29 mm with 60 x 9	0004497800	203,- EUR

Magnetic Stirrers RH basic & digital

IKA®



designed
to work perfectly


IKA® German technology

IKA®-Werke GmbH & Co. KG
Janke & Kunkel-Str. 10 · 79219 Staufen · Deutschland
Tel. +49 7633 831-0 · Fax +49 7633 831-98
sales@ika.de · www.ika.com

RH basic & digital | Highly efficient & economical



Please visit
www.ika.com




RH basic

Adjustable safety temperature

DIN Bushing 12878 for connecting an electronic contact thermometer



Hot Top indicator to prevent burns



RH digital



Digital display for precise monitoring of speed and temperature



DIN Bushing 12878 for connecting an electronic contact thermometer



RH basic/digital with white coated heating plate!

- > Offers excellent chemical resistance
- > Easy to clean

Technical data

Max. stirring quantity (H,O)
Motor rating input / output
Speed range
Speed display
Heat output
Temperature range
Max. magnetic bar (L x Ø)
Safety circuit
Control accuracy with sensor

Heating plate material
Dimensions
Weight

Price

RH basic | RH basic white

15 l
15 / 2 W
50 – 2000 rpm
scale
600 W
50 – 320 °C
80 mm
50 – 370 °C
ETS-D5: ± 0.5 K
ETS-D6: ± 0.2 K

stainless steel 1.4301 | white ceramic
Ø 135 mm
2 kg

350,- EUR | 370,- EUR

Ident. No. 0005019700
Ident. No. 0005029700

RH digital | RH digital white

15 l
15 / 2 W
50 – 2000 rpm
LED
600 W
50 – 320 °C
80 mm
50 – 370 °C
ETS-D5: ± 0.5 K
ETS-D6: ± 0.2 K

stainless steel 1.4301 | white ceramic
Ø 135 mm
2 kg

438,- EUR | 458,- EUR

Ident. No. 0005019800
Ident. No. 0004678000



The newly designed low-cost magnetic stirrers RH basic and digital are now available with enhanced features ensuring better performance and exceptional heating solutions. The strong magnetic field and wide speed range ensures usage for volumes up to 15 liters with ease.



IKA+

The convenient carrying handle and the compact design makes the unit ultra-portable and offers easy transportation

Magnetic Stirrers | Accessories



Electronic contact thermometers

	Ident. No.	Price
1 ETS-D 5	0003378000	281,- EUR
Electronic contact thermometer, -50 – 450 °C, 0.1 K resolution		
2 ETS-D 6	0003378100	361,- EUR
Electronic contact thermometer. Similar to ETS-D5, additionally comes equipped with integrated pH measuring instrument (without pH electrode)		

Temperature sensors for ETS-D5 / D6

	Ident. No.	Price
3 H 62.51	0002735451	132,- EUR
Temperature sensor, stainless steel, Ø 3 mm, 260 mm length		
4 H 66.51	0002735551	195,- EUR
Temperature sensor, stainless steel, glass-coated, Ø 6 mm, 260 mm length		
H 66.53	0004499900	165,- EUR
Temperature sensor, coated with SafeCoat, Ø 3 mm, 260 mm length		

ETS-D5 / D6 Accessories

	Ident. No.	Price
5 H 70	0002735600	128,- EUR
Extension cable, 1 m		

Oil bath attachments

	Ident. No.	Price
H 29	0002829400	123,- EUR
Oil bath attachment, 1 l, aluminum, inner Ø 136-180 mm, 81 mm height for use with oil only		
H 30	0002829500	143,- EUR
Oil bath attachment, 1.5 l, aluminum, inner Ø 136-190 mm, 110 mm height for use with oil only		

Beakers

	Ident. No.	Price
6 H 1000	0004444401	259,- EUR
Beaker, stainless steel, 1 l, inner Ø 160 mm, 79 mm height		
7 H 1500	0004444501	250,- EUR
Beaker, stainless steel, 1.5 l, inner Ø 140 mm, 160 mm height		
8 H 3000	0004444503	251,- EUR
Beaker, stainless steel, 3 l, inner Ø 180 mm, 119 mm height		

Protective covers

	Ident. No.	Price
H 102	0004281600	12,- EUR
Cover for RH basic		
H 103	0004299100	12,- EUR
Cover for RH digital		

Other accessories

	Ident. No.	Price
9 H 16 V	0001545100	45,- EUR
Support rod, Ø 10 mm, 450 mm length		
10 H 16.1	0005000500	12,10 EUR
Extension for support rod		
10 H 44	0002437700	20,- EUR
Boss head clamp		
12 H 38	0003547700	20,- EUR
Holding rod		

